

BONE SCREWS



CAT. No. : BS.05 / TBS.05
Cortex Screws Dia. 3.5mm
 (Stainless Steel & Titanium)
 (Pitch : 1.25mm)

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------

10mm	BS.05.10	TBS.05.10
12mm	BS.05.12	TBS.05.12
14mm	BS.05.14	TBS.05.14
16mm	BS.05.16	TBS.05.16
18mm	BS.05.18	TBS.05.18
20mm	BS.05.20	TBS.05.20
22mm	BS.05.22	TBS.05.22
24mm	BS.05.24	TBS.05.24
26mm	BS.05.26	TBS.05.26
28mm	BS.05.28	TBS.05.28
30mm	BS.05.30	TBS.05.30
32mm	BS.05.32	TBS.05.32
34mm	BS.05.34	TBS.05.34
36mm	BS.05.36	TBS.05.36
38mm	BS.05.38	TBS.05.38
40mm	BS.05.40	TBS.05.40
42mm	BS.05.42	TBS.05.42
44mm	BS.05.44	TBS.05.44
46mm	BS.05.46	TBS.05.46
48mm	BS.05.48	TBS.05.48
50mm	BS.05.50	TBS.05.50



CAT. No. : BS.03 / TBS.03
Cortex Screws Dia. 2.7mm
 (Stainless Steel & Titanium)

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------

6mm	BS.03.06	TBS.03.06
8mm	BS.03.08	TBS.03.08
10mm	BS.03.10	TBS.03.10
12mm	BS.03.12	TBS.03.12
14mm	BS.03.14	TBS.03.14
16mm	BS.03.16	TBS.03.16
18mm	BS.03.18	TBS.03.18
20mm	BS.03.20	TBS.03.20
22mm	BS.03.22	TBS.03.22
24mm	BS.03.24	TBS.03.24
26mm	BS.03.26	TBS.03.26
28mm	BS.03.28	TBS.03.28
30mm	BS.03.30	TBS.03.30

CAT. No. : BS.07 / TBS.07

Cortex Screws Dia. 3.5mm
 (Stainless Steel & Titanium) (Pitch : 1.75mm)

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------

10mm	BS.07.10	TBS.07.10
12mm	BS.07.12	TBS.07.12
14mm	BS.07.14	TBS.07.14
16mm	BS.07.16	TBS.07.16
18mm	BS.07.18	TBS.07.18
20mm	BS.07.20	TBS.07.20
22mm	BS.07.22	TBS.07.22
24mm	BS.07.24	TBS.07.24
26mm	BS.07.26	TBS.07.26
28mm	BS.07.28	TBS.07.28
30mm	BS.07.30	TBS.07.30
32mm	BS.07.32	TBS.07.32
34mm	BS.07.34	TBS.07.34
36mm	BS.07.36	TBS.07.36
38mm	BS.07.38	TBS.07.38
40mm	BS.07.40	TBS.07.40
42mm	BS.07.42	TBS.07.42
44mm	BS.07.44	TBS.07.44
46mm	BS.07.46	TBS.07.46
48mm	BS.07.48	TBS.07.48
50mm	BS.07.50	TBS.07.50



CAT. No. : BS.22 / TBS.22
Washer 10.0mm for Small Screws
 (Stainless Steel & Titanium)
 10.0mm / 4.6mm

OUTER Ø	STAINLESS STEEL	TITANIUM
---------	-----------------	----------

7	BS.22.07	TBS.22.07
13	BS.22.13	TBS.22.13

BONE SCREWS



CAT. No. : BS.09 / TBS.09
4.5mm Cortex Screw
(Pitch : 1.75mm)

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------


14mm	BS.09.14	TBS.09.14
16mm	BS.09.16	TBS.09.16
18mm	BS.09.18	TBS.09.18
20mm	BS.09.20	TBS.09.20
22mm	BS.09.22	TBS.09.22
24mm	BS.09.24	TBS.09.24
26mm	BS.09.26	TBS.09.26
28mm	BS.09.28	TBS.09.28
30mm	BS.09.30	TBS.09.30
32mm	BS.09.32	TBS.09.32
34mm	BS.09.34	TBS.09.34
36mm	BS.09.36	TBS.09.36
38mm	BS.09.38	TBS.09.38
40mm	BS.09.40	TBS.09.40
42mm	BS.09.42	TBS.09.42
44mm	BS.09.44	TBS.09.44
46mm	BS.09.46	TBS.09.46
48mm	BS.09.48	TBS.09.48
50mm	BS.09.50	TBS.09.50
52mm	BS.09.52	TBS.09.52
54mm	BS.09.54	TBS.09.54
56mm	BS.09.56	TBS.09.56
58mm	BS.09.58	TBS.09.58
60mm	BS.09.60	TBS.09.60



CAT. No. : BS.11 / TBS.11
4.0mm Cancellous Screw
(Short Th.)
(Pitch : 1.75mm)

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------

10mm	BS.11.10	TBS.11.10
12mm	BS.11.12	TBS.11.12
14mm	BS.11.14	TBS.11.14
16mm	BS.11.16	TBS.11.16
18mm	BS.11.18	TBS.11.18
20mm	BS.11.20	TBS.11.20
22mm	BS.11.22	TBS.11.22
24mm	BS.11.24	TBS.11.24
26mm	BS.11.26	TBS.11.26
28mm	BS.11.28	TBS.11.28
30mm	BS.11.30	TBS.11.30
35mm	BS.11.35	TBS.11.35
40mm	BS.11.40	TBS.11.40
45mm	BS.11.45	TBS.11.45
50mm	BS.11.50	TBS.11.50
55mm	BS.11.55	TBS.11.55
60mm	BS.11.60	TBS.11.60



CAT. No. : BS.12 / TBS.12
4.0mm Cancellous Screw
(Full Th.)
(Pitch : 1.75mm)

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------

10mm	BS.12.10	TBS.12.10
12mm	BS.12.12	TBS.12.12
14mm	BS.12.14	TBS.12.14
16mm	BS.12.16	TBS.12.16
18mm	BS.12.18	TBS.12.18
20mm	BS.12.20	TBS.12.20
22mm	BS.12.22	TBS.12.22
24mm	BS.12.24	TBS.12.24
26mm	BS.12.26	TBS.12.26
28mm	BS.12.28	TBS.12.28
30mm	BS.12.30	TBS.12.30
35mm	BS.12.35	TBS.12.35
40mm	BS.12.40	TBS.12.40
45mm	BS.12.45	TBS.12.45
50mm	BS.12.50	TBS.12.50
55mm	BS.12.55	TBS.12.55
60mm	BS.12.60	TBS.12.60


BONE SCREWS



CAT. No. : BS.13 / TBS.13
Cancellous Screw Dia. 6.5mm
 (16mm Thread)

Profile

Threaded Diameter : 6.5mm
 Threaded Length : 16.0mm
 Drill Bit for Threading Hole : 3.2mm
 Shaft Diameter : 4.5mm
 Drill Bit for Shaft Optional : 4.5mm
 Core Diameter : 3.0mm
 Hexagonal Socket Width : 3.5mm
 Head Diameter : 8.0mm



CAT. No. : BS.14 / TBS.14
Cancellous Screw
Dia. 6.5mm
 (32mm Thread)

Profile

Threaded Diameter : 6.5mm
 Threaded Length : 32.0mm
 Drill Bit for Threading Hole : 3.2mm
 Shaft Diameter : 4.5mm
 Drill Bit for Shaft Optional : 4.5mm
 Core Diameter : 3.0mm
 Hexagonal Socket Width : 3.5mm
 Head Diameter : 8.0mm



CAT. No. : BS.15 / TBS.15
Cancellous Screw Dia. 6.5mm
 (Full Thread)

Profile

Threaded Diameter : 6.5mm
 Drill Bit for Threading Hole : 3.2mm
 Core Diameter : 3.0mm
 Hexagonal Socket Width : 3.5mm
 Head Diameter : 8.0mm

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------

45mm	BS.14.45	TBS.14.45
50mm	BS.14.50	TBS.14.50
55mm	BS.14.55	TBS.14.55
60mm	BS.14.60	TBS.14.60
65mm	BS.14.65	TBS.14.65
70mm	BS.14.70	TBS.14.70
75mm	BS.14.75	TBS.14.75
80mm	BS.14.80	TBS.14.80
85mm	BS.14.85	TBS.14.85
90mm	BS.14.90	TBS.14.90
95mm	BS.14.95	TBS.14.95
100mm	BS.14.100	TBS.14.100
105mm	BS.14.105	TBS.14.105
110mm	BS.14.110	TBS.14.110
115mm	BS.14.115	TBS.14.115
120mm	BS.14.120	TBS.14.120

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------

20mm	BS.13.20	TBS.13.20
25mm	BS.13.25	TBS.13.25
30mm	BS.13.30	TBS.13.30
35mm	BS.13.35	TBS.13.35
40mm	BS.13.40	TBS.13.40
45mm	BS.13.45	TBS.13.45
50mm	BS.13.50	TBS.13.50
55mm	BS.13.55	TBS.13.55
60mm	BS.13.60	TBS.13.60
65mm	BS.13.65	TBS.13.65
70mm	BS.13.70	TBS.13.70
75mm	BS.13.75	TBS.13.75
80mm	BS.13.80	TBS.13.80
85mm	BS.13.85	TBS.13.85
90mm	BS.13.90	TBS.13.90
95mm	BS.13.95	TBS.13.95
100mm	BS.13.100	TBS.13.100
105mm	BS.13.105	TBS.13.105
110mm	BS.13.110	TBS.13.110
115mm	BS.13.115	TBS.13.115
120mm	BS.13.120	TBS.13.120

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------

20mm	BS.15.20	TBS.15.20
25mm	BS.15.25	TBS.15.25
30mm	BS.15.30	TBS.15.30
35mm	BS.15.35	TBS.15.35
40mm	BS.15.40	TBS.15.40
45mm	BS.15.45	TBS.15.45
50mm	BS.15.50	TBS.15.50
55mm	BS.15.55	TBS.15.55
60mm	BS.15.60	TBS.15.60
65mm	BS.15.65	TBS.15.65
70mm	BS.15.70	TBS.15.70
75mm	BS.15.75	TBS.15.75
80mm	BS.15.80	TBS.15.80
85mm	BS.15.85	TBS.15.85
90mm	BS.15.90	TBS.15.90
95mm	BS.15.95	TBS.15.95
100mm	BS.15.100	TBS.15.100
105mm	BS.15.105	TBS.15.105
110mm	BS.15.110	TBS.15.110
115mm	BS.15.115	TBS.15.115
120mm	BS.15.120	TBS.15.120

BONE SCREWS



CAT. No. : BS.23 / TBS.23
Malleolar Screws Dia. 4.5mm

Profile

Threaded Diameter : 4.5mm
Drill Bit for Threading Hole : 3.2mm
Shaft Diameter : 3.0mm
Core Diameter : 3.0mm
Hexagonal Socket Width : 3.5mm
Head Diameter : 8.0mm

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------

20mm	BS.23.20	TBS.23.20
25mm	BS.23.25	TBS.23.25
30mm	BS.23.30	TBS.23.30
35mm	BS.23.35	TBS.23.35
40mm	BS.23.40	TBS.23.40
45mm	BS.23.45	TBS.23.45
50mm	BS.23.50	TBS.23.50
55mm	BS.23.55	TBS.23.55
60mm	BS.23.60	TBS.23.60
65mm	BS.23.65	TBS.23.65
70mm	BS.23.70	TBS.23.70
75mm	BS.23.75	TBS.23.75
80mm	BS.23.80	TBS.23.80
85mm	BS.23.85	TBS.23.85
90mm	BS.23.90	TBS.23.90



CAT. No. : BS.16 / TBS.16
Cannulated Cancellous Screw
Dia. 4.0mm
(Short Thread)

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------

10mm	BS.16.10	TBS.16.10
12mm	BS.16.12	TBS.16.12
14mm	BS.16.14	TBS.16.14
16mm	BS.16.16	TBS.16.16
18mm	BS.16.18	TBS.16.18
20mm	BS.16.20	TBS.16.20
22mm	BS.16.22	TBS.16.22
24mm	BS.16.24	TBS.16.24
26mm	BS.16.26	TBS.16.26
28mm	BS.16.28	TBS.16.28
30mm	BS.16.30	TBS.16.30
35mm	BS.16.35	TBS.16.35
40mm	BS.16.40	TBS.16.40
45mm	BS.16.45	TBS.16.45
50mm	BS.16.50	TBS.16.50
55mm	BS.16.55	TBS.16.55
60mm	BS.16.60	TBS.16.60



CAT. No. : BS.17 / TBS.17
Cannulated Cancellous Screw
Dia. 4.0mm
(Fully Thread)

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------

10mm	BS.17.10	TBS.17.10
12mm	BS.17.12	TBS.17.12
14mm	BS.17.14	TBS.17.14
16mm	BS.17.16	TBS.17.16
18mm	BS.17.18	TBS.17.18
20mm	BS.17.20	TBS.17.20
22mm	BS.17.22	TBS.17.22
24mm	BS.17.24	TBS.17.24
26mm	BS.17.26	TBS.17.26
28mm	BS.17.28	TBS.17.28
30mm	BS.17.30	TBS.17.30

LENGTH	STAINLESS STEEL	TITANIUM
--------	-----------------	----------

35mm	BS.17.35	TBS.17.35
40mm	BS.17.40	TBS.17.40
45mm	BS.17.45	TBS.17.45
50mm	BS.17.50	TBS.17.50
55mm	BS.17.55	TBS.17.55
60mm	BS.17.60	TBS.17.60

BONE SCREWS



CAT. No. : BS.18 / TBS.18
Cannulated Cancellous Screw
Dia. 6.5mm
 (16mm Thread)

LENGTH	STAINLESS STEEL	TITANIUM
30mm	BS.18.30	TBS.18.30
35mm	BS.18.35	TBS.18.35
40mm	BS.18.40	TBS.18.40
45mm	BS.18.45	TBS.18.45
50mm	BS.18.50	TBS.18.50
55mm	BS.18.55	TBS.18.55
60mm	BS.18.60	TBS.18.60
65mm	BS.18.65	TBS.18.65
70mm	BS.18.70	TBS.18.70
75mm	BS.18.75	TBS.18.75
80mm	BS.18.80	TBS.18.80
85mm	BS.18.85	TBS.18.85
90mm	BS.18.90	TBS.18.90
95mm	BS.18.95	TBS.18.95
100mm	BS.18.100	TBS.18.100
105mm	BS.18.105	TBS.18.105
110mm	BS.18.110	TBS.18.110



CAT. No. : BS.19 / TBS.19
Cannulated Cancellous Screw
Dia. 6.5mm
 (32mm Thread)

LENGTH	STAINLESS STEEL	TITANIUM
40mm	BS.18.40	TBS.18.40
45mm	BS.18.45	TBS.18.45
50mm	BS.18.50	TBS.18.50
55mm	BS.18.55	TBS.18.55
60mm	BS.18.60	TBS.18.60
65mm	BS.18.65	TBS.18.65
70mm	BS.18.70	TBS.18.70
75mm	BS.18.75	TBS.18.75
80mm	BS.18.80	TBS.18.80
85mm	BS.18.85	TBS.18.85
90mm	BS.18.90	TBS.18.90
95mm	BS.18.95	TBS.18.95
100mm	BS.18.100	TBS.18.100
105mm	BS.18.105	TBS.18.105
110mm	BS.18.110	TBS.18.110