

HIP SYSTEM



CAT. No. : HS.04
Bipolar Prosthesis
Non-Fenestrated (Sterile)

LENGTH	STAINLESS STEEL	HEAD Ø
155mm	HS.04.37	37.0mm
155mm	HS.04.39	39.0mm
155mm	HS.04.41	41.0mm
155mm	HS.04.43	43.0mm
155mm	HS.04.45	45.0mm
155mm	HS.04.47	47.0mm
155mm	HS.04.49	49.0mm
155mm	HS.04.51	51.0mm
155mm	HS.04.53	53.0mm
155mm	HS.04.55	55.0mm

LENGTH	STAINLESS STEEL	HEAD Ø
140mm	HS.01.37	37.0mm
140mm	HS.01.39	39.0mm
140mm	HS.01.41	41.0mm
140mm	HS.01.43	43.0mm
140mm	HS.01.45	45.0mm
140mm	HS.01.47	47.0mm
140mm	HS.01.49	49.0mm
140mm	HS.01.51	51.0mm
140mm	HS.01.53	53.0mm
140mm	HS.01.55	55.0mm

CAT. No. : HS.01
Austin Moore
Endo-Prosthesis
(Sterile)



HIP SYSTEM



CAT. No. : BP.42
Dynamic Hip Plate 125°

HOLE	STAINLESS STEEL	LENGTH
2 Hole	BP.42.02	46mm
4 Hole	BP.42.04	78mm
5 Hole	BP.42.05	94mm
6 Hole	BP.42.06	110mm
8 Hole	BP.42.08	142mm
10 Hole	BP.42.10	174mm

CAT. No. : BP.43
Dynamic Hip Plate 130°

HOLE	STAINLESS STEEL	LENGTH
2 Hole	BP.43.02	46mm
4 Hole	BP.43.04	78mm
5 Hole	BP.43.05	94mm
6 Hole	BP.43.06	110mm
8 Hole	BP.43.08	142mm
10 Hole	BP.43.10	174mm
12 Hole	BP.43.12	206mm
14 Hole	BP.43.14	238mm
16 Hole	BP.43.16	270mm

CAT. No. : BP.44
Dynamic Hip Plate 135°

HOLE	STAINLESS STEEL	LENGTH
2 Hole	BP.44.02	46mm
4 Hole	BP.44.04	78mm
5 Hole	BP.44.05	94mm
6 Hole	BP.44.06	110mm
8 Hole	BP.44.08	142mm
10 Hole	BP.44.10	174mm
12 Hole	BP.44.12	206mm
14 Hole	BP.44.14	238mm
16 Hole	BP.44.16	270mm

CAT. No. : BP.45
Dynamic Hip Plate 140°


HOLE	STAINLESS STEEL	LENGTH
2 Hole	BP.45.02	46mm
4 Hole	BP.45.04	78mm
5 Hole	BP.45.05	94mm
6 Hole	BP.45.06	110mm
8 Hole	BP.45.08	142mm
10 Hole	BP.45.10	174mm
12 Hole	BP.45.12	206mm
14 Hole	BP.45.14	238mm
16 Hole	BP.45.16	270mm

CAT. No. : BP.46
Dynamic Hip Plate 145°

HOLE	STAINLESS STEEL	LENGTH
2 Hole	BP.46.02	46mm
4 Hole	BP.46.04	78mm
5 Hole	BP.46.05	94mm
6 Hole	BP.46.06	110mm
8 Hole	BP.46.08	142mm
10 Hole	BP.46.10	174mm
12 Hole	BP.46.12	206mm
14 Hole	BP.46.14	238mm
16 Hole	BP.46.16	270mm

CAT. No. : BP.47
Dynamic Hip Plate 150°

HOLE	STAINLESS STEEL	LENGTH
2 Hole	BP.50.02	46mm
4 Hole	BP.50.04	78mm
5 Hole	BP.50.05	94mm
6 Hole	BP.50.06	110mm
8 Hole	BP.50.08	142mm
10 Hole	BP.50.10	174mm
12 Hole	BP.50.12	206mm
14 Hole	BP.50.14	238mm
16 Hole	BP.50.16	270mm



CAT. No. : BP.40
DCS Plate DCP Hole 95°

HOLE	STAINLESS STEEL	LENGTH
4 Hole	BP.40.04	81mm
5 Hole	BP.40.05	97mm
6 Hole	BP.40.06	113mm
8 Hole	BP.40.08	145mm
10 Hole	BP.40.10	177mm
12 Hole	BP.40.12	209mm
14 Hole	BP.40.14	241mm
16 Hole	BP.40.16	273mm
18 Hole	BP.40.18	305mm
20 Hole	BP.40.20	337mm
22 Hole	BP.40.22	369mm



CAT. No. : CS.29
DHS / DCS Compression Screw

LENGTH	STAINLESS STEEL
36mm	CS.29.36



CAT. No. : BS.21
DHS / DCS Screw 12.5mm

LENGTH	STAINLESS STEEL
50mm	BS.21.50
55mm	BS.21.55
60mm	BS.21.60
65mm	BS.21.65
70mm	BS.21.70
75mm	BS.21.75
80mm	BS.21.80
85mm	BS.21.85
90mm	BS.21.90
95mm	BS.21.95
100mm	BS.21.100
105mm	BS.21.105
110mm	BS.21.110

